



DEVICE SPECIFICATIONS




Binocular, OLED color display that integrates prescription lenses.



21.5 MP camera with high acuity auto-focus liquid lens, and 24x zoom.



Accelerometers, gyroscope, IR range finder and temperature sensor.



Advanced UI with HDMI, Bluetooth®, Wi-Fi, and removable SD cards.

DISPLAY SPECIFICATIONS	
TYPE	Dual independent high contrast OLEDs.
RESOLUTION	1024 X 768
INTER-PUPIL DISTANCE	Adjustable
FIELD OF VIEW	Diagonal 37.5°
FRAME ATTACHMENT	Magnetic (patent pending)

HMD SPECIFICATIONS	
TYPE	Bi-optic (patented)
ADJUSTABLE BIOPTIC TILT	+45° to -10°
MASS	104g
SYSTEM	Android (5.0)
SYSTEM ACUITY	>20/20

