


## The best enhanced vision device is ready for life on the go



### Device Specifications

- Streamlined visor with adjustable arms fits over prescription glasses
- Patented bioptic tilt adjusts between +35° to -10°
- Adjustable halo band securely fits children and adults
- 2 rechargeable lithium-ion batteries, each with up to 3 hours of continuous use
- 256 GB storage
- 3 integrated speakers
- Built-in flashlight

### Display Specifications

- System acuity > 20/20
- 2 dual independent 1280X960 high resolution OLED color screens
- 18 MP camera for stunning clarity
- Advanced software auto-focus, image stabilization, and personalization options
- Adjustable inter-pupil distance
- Diagonal 37.5° field of view

### Connectivity and Controls

- Touchpad and menu control integrated into headset
- Wireless remote and Web account controls
- Android® and Apple® apps
- Wi-Fi, Bluetooth and HDMI connectivity
- eCast and eMirror with the eSight mobile apps
- eSupport and Shared Accounts

See new possibilities, with eSight 4.